

TRENGUARD HYBRID -POTILASTUKI TRENDELENBURG-ASENTOON

TrenGuard on edistyksellinen ratkaisu potilaan asetteluun MIS- ja robottiaavusteisessa kirurgiassa.

Hyödyt:


- Ehkäisee potilaan liukumista Trendelenburg-asennossa erityisen niskakuoppapehmusteen avulla.
- Ei tarvetta olkatuille. Pehmeät lateraaliset stabilointityynyntä kontrolloivat vartalon painopisteen muutoksia.
- Ei tarvetta vartalohihnoille helpottaen potilaan ventilaatiota.
- Tehokas ja helppokäyttöinen. Sopii kaikkiin Merivaaran leikkauspöytiin (myydään vain leikkauspöydän mukana).
- Vähemmän jätettä uudelleen käytettävien lateraalisten stabilointityynyjen ansiosta.
- Alhaisemmat materiaalikustannukset.
- Järjestelmä sopii jopa 250 kg painoisille potilaille.





TrenGuard 450 Hybrid starter package Standard


- Including:
- Patient Support Frame
 - Clamps for EUR rails
 - Single-use soft head stabilizing pillow and Speed Bump bolster
 - Lateral Stabilizing Pillows



TrenGuard single-use wedge package

For bariatric patients and those without discernible cervical notch.

- Including:
- 12 single-use integrated wedge and head stabilizing pillows
 - Order No. 56401



TrenGuard single-use package

- Including:
- 12 single-use Speed Bump bolsters which prevent patient sliding
 - 12 single-use soft head stabilizing pillows which control side-to-side rocking of the head
 - Order No. 55201



TrenGuard single-use package

For bariatric patients.

- Including:
- 12 single-use Speed Bump bolsters which prevent patient sliding
 - 12 single-use soft head stabilizing pillows which control side-to-side rocking of the head
 - Order No. 56201



TrenGuard reusable package

- Including:
- One reusable soft lateral stabilizing pillow control body mass shift during transition to Trendelenburg
 - Order No. 55251



TrenGuard reusable package

For bariatric patients.

- Including:
- One reusable soft lateral stabilizing pillow to control body mass shift during transition to Trendelenburg
 - Order No. 56251