

# DORO HORSESHOE SYSTEM

Universal Horseshoe Headrest for adult or pediatric applications in prone or supine positioning. Mounted to a Doro Swivel Adaptor Aluminium or to the Doro Horseshoe Table Adaptor. Gel Pads interchangeable for adult or pediatric application. Permits vertical and lateral adjustments for proper patient positioning.



## Horseshoe system

100 300 900 Adult headrest with Extension Bar  
 100 300 800 Pediatric headrest with Extension Bar

## Pads

100300910 Adult Right  
 100300911 Adult Left  
 100300810 Pediatric Right  
 100300811 Pediatric Left

## Cervical Spine Support

100301200

## Tuotekoodi

100 300 900 100 300 800 100300910  
 100300911 100300810 100300811  
 100301200

---